

IPPK antibody



Catalog Number: 111821

Product name

IPPK antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

IPPK Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IHC detected in human gliomas. Observed molecular weight by Western-blot: 56kd

Preparation

This antibody was obtained by immunization of IPPK recombinant protein (Accession Number: NM_022755). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111821(IPPK antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:111821(IPPK antibody) at dilution of 1:100 (under 10x lens)