

# Anti-KCNH5 antibody

Catalog Number: 175821

## Product name

Anti-KCNH5 antibody

WB: 1:1000-1:2000

IHC-P: 1:100-200

## Specificity

Human, Rat, Mouse

## Antibody description

Rabbit polyclonal antibody to KCNH5

## Preparation

This antigen of this antibody was synthetic peptide of kv10.2

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P

## Dilutions

## Validations




Fig1::; Western blot analysis of 1) Rat Brain Tissue-Low Molecular Protein Marker, 2)Rat Brain Tissue-High Molecular Protein Marker, 3) Mouse Brain Tissue-Low Molecular Protein Marker, 4) Mouse Brain Tissue- High Molecular Protein Marker, 5) HepG2-Low Molecular P




Fig2::; Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using Kv10.2 Rabbit pAb diluted at 1:200.