

Lexatumumab reference antibody



Catalog Number: 176821

Product name

Lexatumumab reference antibody

Specificity

Human TNFRSF10B / TRAILR2 / CD262

Antibody description

Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (Lexatumumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (lexatumumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (lexatumumab) is 99.04%, determined by



SEC-HPLC.

Bioactivity: ELISA

Immobilized human TNFRSF10B/TRAILR2 / CD262, Fc Tag at 2 µg/mL can bind lexatumumab, EC50=0.03685 µg/mL.