

GPR183/EBI2 antibody



Catalog Number: 111121

Product name

GPR183/EBI2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GPR183/EBI2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, Jurkat cells, mouse lung tissue, Raji cells. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 55-60 kDa

Preparation

This antibody was obtained by immunization of GPR183/EBI2 recombinant protein (Accession Number: XM_017020405). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111121(GPR183 antibody) at dilution of 1:300



IP Result of anti-GPR183 (IP:Catalog No:111121, 4ug; Detection:Catalog No:111121 1:500) with mouse lung tissue lysate 4000ug.