

NEURL antibody

Catalog Number: 113121

Product name

NEURL antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

NEURL Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, MCF7 cells, Y79 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 60-62 kDa

Preparation

This antibody was obtained by immunization of NEURL recombinant protein (Accession Number: NM_004210). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:113121 (NEURL antibody) at dilution of 1:800



IP Result of anti-NEURL (IP:Catalog No:113121, 4ug; Detection:Catalog No:113121 1:600) with mouse brain tissue lysate 2560ug.