

# TID1 antibody

Catalog Number: 116121

## Product name

TID1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TID1 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human pancreas cancer tissue. Positive IP detected in K-562 cells. Positive WB detected in Jurkat cells, K-562 cells. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of TID1 recombinant protein (Accession Number: NM\_005147). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116121(DNAJA3 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:116121(DNAJA3 antibody) at dilution of 1:50



(under 40x lens)

IP Result of anti-DNAJA3 (IP:Catalog No:116121, 4ug; Detection:Catalog No:116121 1:500) with

# TID1 antibody



Catalog Number: 116121

---

K-562 cells lysate 3200ug.