

Anti-IMP3 (1F1) Mouse antibody



Catalog Number: 168121

Product name

Anti-IMP3 (1F1) Mouse antibody

Specificity

Human

Antibody description

IMP3 (1F1) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human IMP3 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of IMP3 in MCF7,SKVO3,NTERA2,Jurkat,Hela and K562 cell lysates using IMP3 mouse mAb (1:1000 diluted).Predicted band size:70KDa.Observed band size:70KDa.