

AP3B1 antibody

Catalog Number: 108121

Product name

AP3B1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP3B1 Rabbit Polyclonal antibody. Positive IP detected in COLO 320 cells. Positive WB detected in mouse thymus tissue, COLO 320 cells. Observed molecular weight by Western-blot: 140 kDa

Preparation

This antibody was obtained by immunization of AP3B1 recombinant protein (Accession Number: NM_003664). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:108121(AP3B1 antibody) at dilution of 1:500



IP Result of anti-AP3B1 (IP:Catalog No:108121, 3ug; Detection:Catalog No:108121 1:500) with COLO 320 cells lysate 2500ug.