

Anti-CD55 Mouse antibody

Catalog Number: 169121

Product name

Anti-CD55 Mouse antibody

Specificity

Human

Antibody description

CD55 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CD55 (AA: 238-357) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

FC


Validations



Flow cytometric analysis of K562 cells using CD55 mouse mAb (green) and negative control (red).



Western blot analysis using CD55 mAb against human CD55 recombinant protein. (Expected MW



is 38.5 kDa)

Western blot analysis using CD55 mouse mAb against Raji (1) and K562 (2) cell lysate.