

# FANCL antibody



Catalog Number: 110521

---

## Product name

FANCL antibody

## Immunogen

[Human FANCL Recombinant protein \(His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FANCL Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells. Observed molecular weight by Western-blot: 41-43kd

## Preparation

This antibody was obtained by immunization of FANCL recombinant protein (Accession Number: NM\_018062). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110521(FANCL antibody) at dilution of 1:300