

Anti-SPN antibody



Catalog Number: 142521

Product name

Anti-SPN antibody

Specificity

Human SPN / CD43

Antibody description

Rabbit Polyclonal to Human SPN

Preparation

Produced in rabbits immunized with a synthetic peptide corresponding to the N-terminus of the Human SPN / CD43, and purified by antigen affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It

is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IP

Dilutions

IP: 4-8 µl/mg of lysate

Validations



Items	Lane	A
Sample (whole cell lysate)		Jurkat
Sample quantity		0.5 mg
IP antibody quantity		4 µl
Protein G agarose		15 µl of 50% Protein G agarose
Gel		13% SDS-PAGE reducing gel
Primary antibody		HBcAG-2-SMCC-SINO1324 antibody at 10 µg/ml
Secondary antibody		Clean-Blot IP Detection Reagent (HRP) at 1:500 dilution

Human SPN/CD43 Immunoprecipitation(IP) 14617