

Anti-E Cadherin Rabbit antibody



Catalog Number: 166921

Product name

Anti-E Cadherin Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Rat

Applications

WB, IP

Antibody description

E Cadherin Rabbit polyclonal antibody

Dilutions

WB: 1/1000

IP: 1/20

Preparation

Antigen: Recombinant protein of mouse E Cadherin

Validations

Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal



Western blot detection of E Cadherin in Mouse brain lysates using E Cadherin Rabbit pAb(1:1000 diluted). Predicted band size: 98(cleavages), 135KDa. Observed band size: 80-120(cleavages), 135KDa.