

# Anti-KAT7/Hbo1/MYST2 Rabbit antibody



Catalog Number: 167921

## Product name

Anti-KAT7/Hbo1/MYST2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

KAT7/Hbo1/MYST2 Rabbit monoclonal antibody

## Preparation

Antigen: Recombinant protein of human KAT7/Hbo1/MYST2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, IP

## Dilutions


WB: 1/2000

IHC: 1/50

ICC/IF: 1/200

IP: 1/20

## Validations



Western blot detection of KAT7/Hbo1/MYST2 in HeLa, A549, MCF-7, 293T cell lysates using KAT7/Hbo1/MYST2 Rabbit mAb(1:1000 diluted). Predicted band size: 71 kDa. Observed band size: 71 kDa.