

Anti-CRTC3 Rabbit antibody



Catalog Number: 168921

Product name

Anti-CRTC3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

CRTC3 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human CRTC3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/2000

ICC/IF: 1/20

Validations



Western blot detection of CRTC3 in K562,Rat Brain,C6,3T3,Hela cell lysates using CRTC3 Rabbit mAb(1:1000 diluted).Predicted band size:67kDa.Observed band size:78kDa.