

PROSC antibody



Catalog Number: 114221

Product name

PROSC antibody

Immunogen

[Human PROSC Recombinant protein \(GST tag\)](#)

Specificity

Human; other species not tested.

Antibody description

PROSC Rabbit Polyclonal antibody. Positive IF detected in L02 cells. Positive WB detected in L02 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of PROSC recombinant protein (Accession Number: NM_007198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:114221(PROSC Antibody) at dilution of 1:1000



Immunofluorescent analysis of L02 cells using Catalog No:114221(PROSC Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)