

BCL2L15 Polyclonal Antibody



Catalog Number: 165221

Product name

BCL2L15 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to BCL2L15

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human BCL2L15 (NP_001010922.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:100

Validations



Immunofluorescence - BCL2L15 Polyclonal Antibody

Immunofluorescence analysis of A549 cells using BCL2L15 antibody . Blue: DAPI for nuclear staining.