

# HRP-conjugated Flag Tag antibody



Catalog Number: 107221

## Product name

HRP-conjugated Flag Tag antibody

## Specificity

Human; other species not tested.

## Antibody description

HRP-conjugated Flag Tag Mouse Monoclonal antibody. Positive WB detected in Recombinant protein. Positive IP detected in Recombinant protein. Positive IF detected in Transfected HEK-293 cells.

## Preparation

This antibody was obtained by immunization of Peptide. Purification method: Protein A purified.

## Formulation

PBS with 0.02% Thimerosal, 50% Glycerol, pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b

## Applications

ELISA, WB, IP, IF


## Dilutions

Recommended Dilution:

WB: 1:5000-1:50000

IP: N/AIF: 1:20-1:200

## Validations



Flag tagged fusion protein (20ng/lane) was subjected to SDS-PAGE followed by western blot with HRP-66008 at dilution of 1:40000.