

Mitofusin 2 Polyclonal Antibody



Catalog Number: 164621

Product name

Mitofusin 2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Mitofusin 2

Preparation

Antigen: Recombinant protein of human Mitofusin 2

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IP


Dilutions

WB 1:1000 - 1:4000

IHC 1:50 - 1:200

IP 1:20 - 1:50

Validations



Western blot - Mitofusin 2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using Mitofusin 2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.