

Anti-BRAT1 Rabbit antibody



Catalog Number: 168621

Product name

Rabbit IgG

Anti-BRAT1 Rabbit antibody

Applications

Specificity

WB, FC, IP

Human

Dilutions

Antibody description

WB: 1/1000

BRAT1 Rabbit monoclonal antibody

FC: 1/20

Preparation

IP: 1/20

Antigen: Recombinant protein of human BRAT1

Validations

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA.


Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type



Western blot detection of BRAT1 in HCT116, HeLa cell lysates using BRAT1 Rabbit mAb(1:1000 diluted). Predicted band size: 88kDa. Observed band size: 88kDa.