

LRP2BP antibody



Catalog Number: 112321

Product name

LRP2BP antibody

Immunogen

[Human LRP2BP Recombinant protein \(His tag\)](#)

Specificity

Human, Rat; other species not tested.

Antibody description

LRP2BP Rabbit Polyclonal antibody. Positive WB detected in A549 cells, Jurkat cells, rat lung tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of LRP2BP recombinant protein (Accession Number: NM_001385601). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:112321(LRP2BP Antibody) at dilution of 1:1000