

PTPN4 antibody

Catalog Number: 114321

Product name

PTPN4 antibody

Immunogen

[Human PTPN4 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PTPN4 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, HepG2 cells, human lung tissue, mouse thymus tissue, PC-3 cells. Positive IHC detected in human pancreas cancer tissue. Observed molecular weight by Western-blot: 50-60 kDa

Preparation

This antibody was obtained by immunization of PTPN4 recombinant protein (Accession Number: NM_002830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:114321(PTPN4 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:114321(PTPN4 antibody) at dilution of 1:50 (under 10x lens)