

Anti-HTR6 antibody



Catalog Number: 176321

Product name

Anti-HTR6 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to HTR6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human 5-HT₆. aa range:221-270

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of A549 cells using SR-6 Polyclonal Antibody