

TRIM6 antibody

Catalog Number: 116321

Product name

TRIM6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TRIM6 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in COLO 320 cells, human liver tissue. Observed molecular weight by Western-blot: 56 kDa

Preparation

This antibody was obtained by immunization of TRIM6 recombinant protein (Accession Number: BC065575). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:116321 (TRIM6 antibody) at dilution of 1:50



(under 10x lens)

COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:116321 (TRIM6 Antibody) at dilution of 1:1000