

Human IgG antibody



Catalog Number: 111721

Product name

Human IgG antibody

Specificity

Human; other species not tested.

Antibody description

Human IgG Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 25 kDa, 50 kDa, 150 kDa

Preparation

This antibody was obtained by immunization of Recombinant Protein. Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:2000-1:20000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:111721(Human-IgG Antibody) at dilution of 1:5000