

Product name

ACADS antibody

Specificity

Human, Mouse, Rat, Monkey; other species not tested.

Antibody description

ACADS Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human heart tissue. Positive IP detected in HEK-293 cells. Positive WB detected in human brain tissue, COS-7 cells, HEK-293 cells, HepG2 cells, mouse kidney tissue, mouse liver tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of ACADS recombinant protein (Accession Number: NM_000017). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:200-1:2000

IHC: 1:20-1:200

Validations





human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107721(ACADS antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107721(ACADS antibody) at dilution of 1:100 (under 10x lens)

ACADS antibody

Catalog Number: 107721



Immunohistochemical of paraffin-embedded human kidney using Catalog No:107721(ACADS antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-ACADS (IP:Catalog No:107721, 4ug; Detection:Catalog No:107721 1:500) with HEK-293 cells lysate 3000ug.