

# Anti-AMPK gamma 1 Rabbit antibody



Catalog Number: 168721

## Product name

Rabbit IgG

Anti-AMPK gamma 1 Rabbit antibody

## Applications

## Specificity

WB, FC, IP

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

AMPK gamma 1 Rabbit monoclonal antibody

FC: 1/20

## Preparation

IP: 1/20

Antigen: A synthetic peptide of human AMPK gamma 1

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of AMPK gamma 1 in CHO-K1, C6, Jurkat cell lysates using AMPK gamma 1 Rabbit mAb(1:500 diluted). Predicted band size: 38kDa. Observed band size: 38kDa.