

Anti-hnRNP G Rabbit antibody



Catalog Number: 167031

Product name

Anti-hnRNP G Rabbit antibody

WB: 1/1000

IHC: 1/20-1/100

Specificity

Human, Mouse, Rat

ICC/IF: 1/20-1/50

Antibody description

hnRNP G Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human hnRNP G

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF

Dilutions

Validations



Western blot detection of hnRNP G in K562, HeLa cell lysates using hnRNP G Rabbit pAb(1:1000 diluted). Predicted band size:42kDa. Observed band size:42kDa.



Immunohistochemistry of hnRNP G in paraffin-embedded Human lung cancer tissue using hnRNP G Rabbit pAb at dilution 1/20