

CAPG antibody



Catalog Number: 107031

Product name

CAPG antibody

Specificity

Human; other species not tested.

Antibody description

CAPG Mouse Monoclonal antibody. Positive WB detected in Jurkat cells, HeLa cells, HL-60 cells. Observed molecular weight by Western-blot: 39-40 kDa

Preparation

This antibody was obtained by immunization of CAPG recombinant protein (Accession Number: XM_011533122). Purification method: Protein G purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG1

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:107031(CAPG Antibody) at dilution of 1:1000