

# ZNF750 antibody



Catalog Number: 117031

## Product name

ZNF750 antibody

## Specificity

human, mouse; other species not tested.

## Antibody description

ZNF750 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 70-77kd

## Preparation

This antibody was obtained by immunization of ZNF750 recombinant protein (Accession Number: NM\_024702). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:117031(ZNF750 antibody) at dilution of 1:500