

ZNF750 antibody



Catalog Number: 117031

Product name

ZNF750 antibody

Specificity

human, mouse; other species not tested.

Antibody description

ZNF750 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue. Observed molecular weight by Western-blot: 70-77kd

Preparation

This antibody was obtained by immunization of ZNF750 recombinant protein (Accession Number: NM_024702). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:117031(ZNF750 antibody) at dilution of 1:500