

# F7 antibody

Catalog Number: 110431

---

## Product name

F7 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

F7 Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human plasma tissue, mouse liver tissue. Observed molecular weight by Western-blot: 49-55 kDa

## Preparation

This antibody was obtained by immunization of F7 recombinant protein (Accession Number: NM\_019616). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110431(F7 Antibody) at dilution of 1:300