

Mouse CD99 (Fc Tag) recombinant protein



Catalog Number: 503431

General Information

Gene Name Synonym

MIC2-like protein 1

Protein Construction

A DNA sequence encoding the mouse CD99 (NP_612182.1) extracellular domain (Met 1-Gly 137) was fused with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Thr 21

Molecular Mass

The recombinant mouse CD99/Fc chimera is a

disulfide-linked homodimeric protein. The reduced monomer consists of 350 amino acids and predicts a molecular mass of 37.5 kDa. As a result of glycosylation, the apparent molecular mass of rm CD99/Fc monomer is approximately 47-55 kDa in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse CD99 / PILR1 Protein (Fc Tag) SDS-PAGE