

Anti-CD163L1 antibody

Catalog Number: 176431

Product name

Anti-CD163L1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CD163L1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human cd163l1. aa range:391-440

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of HeLa cells using CD163b Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2.; Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100