

IRF1 antibody



Catalog Number: 111831

Product name

IRF1 antibody

ELISA, WB, IHC

Dilutions

Recommended Dilution:

Specificity

Human, Mouse, Rat; other species not tested.

WB: 1:200-1:1000

IHC: 1:20-1:200


Antibody description

IRF1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, Jurkat cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 45-50 kDa

Validations

Preparation

This antibody was obtained by immunization of IRF1 recombinant protein (Accession Number: NM_002198). Purification method: Antigen affinity purified.



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:111831(IRF1 Antibody) at dilution of 1:300

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:111831(IRF1 Antibody) at dilution of 1:50 (under 40x lens)