

SULT2A1 antibody



Catalog Number: 115831

Product name

SULT2A1 antibody

Specificity

Human; other species not tested.

Antibody description

SULT2A1 Rabbit Polyclonal antibody. Positive WB detected in human adrenal gland tissue, human liver tissue, L02 cells. Observed molecular weight by Western-blot: 34 kDa

Preparation

This antibody was obtained by immunization of SULT2A1 recombinant protein (Accession Number: NM_003167). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human adrenal gland tissue were subjected to SDS PAGE followed by western blot with Catalog No:115831(SULT2A1 antibody) at dilution of 1:1000