

Anti-MVK Rabbit antibody



Catalog Number: 167831

Product name

Anti-MVK Rabbit antibody

Specificity

Human

Antibody description

MVK Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human MVK

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

WB, FC, IP


Dilutions

WB: 1/1000

FC: 1/20

IP: 1/20

Validations



Western blot detection of MVK in HeLa,A549,HL-60,U251,U87-MG cell lysates using MVK Rabbit mAb(1:1000 diluted).Predicted band size:43kDa.Observed band size:43kDa.