

PPIL6 antibody



Catalog Number: 114131

Product name

PPIL6 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

PPIL6 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells. Observed molecular weight by Western-blot: 30kd

Preparation

This antibody was obtained by immunization of PPIL6 recombinant protein (Accession Number: NM_173672). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:114131(PPIL6 antibody) at dilution of 1:500