

Product name

APC antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

APC Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human colon cancer tissue, human colon tissue, human endometrial cancer tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000038).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:108131 (APC Antibody) at dilution of 1:50. Heat mediated antigen retrieved with Citric acid buffer, pH6.0.



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:108131 (APC Antibody) at dilution of 1:50. Heat mediated antigen retrieved with Citric acid buffer, pH6.0.