

# Anti-TLR4 (Center) Rabbit antibody



Catalog Number: 169131

## Product name

Anti-TLR4 (Center) Rabbit antibody

## Specificity

Human

## Antibody description

TLR4 (Center) Rabbit polyclonal antibody

## Preparation

Antigen: This TLR4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 669-698 amino acids from the Central region of human TLR4.

## Formulation

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB, FC


## Dilutions

WB::1:1000FCM::1:10~50

## Validations



Western blot analysis of TLR4 Antibody (Center) (Cat. #169131) in HL-60 cell line lysates (35ug/lane). TLR4 (arrow) was detected using the



purified Pab.

TLR4 Antibody (Center) (Cat. #169131) flow cytometric analysis of U251 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis.