

# HOXB1 antibody



Catalog Number: 111531

## Product name

HOXB1 antibody

## Immunogen

[Human HOXB1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HOXB1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 41-44 kDa

## Preparation

This antibody was obtained by immunization of HOXB1 recombinant protein (Accession Number: BC096192). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111531(HOXB1 Antibody) at dilution of 1:300