

# Pan Acetyl-Lysine mAb



Catalog Number: 162531

---

## Product name

Pan Acetyl-Lysine mAb

## Specificity

ALL

## Antibody description

Monoclonal antibody to Pan Acetyl-Lysine

## Preparation

Antigen: Acetylated BSA

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB

## Dilutions

WB 1:500 - 1:1000

## Validations



Western blot - Pan Acetyl-Lysine mAb

Western blot analysis of extracts of HeLa cells, using Pan Acetyllysine antibody .