

Human PDRG1 / C20orf126 (His Tag) recombinant protein



Catalog Number: 502531

General Information

Protein Construction

A DNA sequence encoding the human PDRG1 (Q9NUG6) (Met1-Gly133) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human PDRG1 consists of 148 amino acids and predicts a molecular mass of 17.3 KDa. It migrates as an approximately 12 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, 10% Glycerol, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human PDRG1 / C20orf126 Protein (His Tag)
SDS-PAGE