

SPAG4L antibody



Catalog Number: 115531

Product name

SPAG4L antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

SPAG4L Rabbit Polyclonal antibody. Positive WB detected in human testis tissue, HeLa cells. Observed molecular weight by Western-blot: 43 kDa

Preparation

This antibody was obtained by immunization of SPAG4L recombinant protein (Accession Number: NM_080675). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:115531(SPAG4L antibody) at dilution of 1:400