

UBE2L6 antibody



Catalog Number: 116531

Product name

UBE2L6 antibody

IP: 1:500-1:5000

IF: 1:10-1:100

Specificity

Human; other species not tested.

Antibody description

UBE2L6 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Positive IP detected in Jurkat cells. Positive IF detected in Hela cells. Observed molecular weight by Western-blot: 18 kDa

Preparation

This antibody was obtained by immunization of UBE2L6 recombinant protein (Accession Number: NM_004223). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116531(UBE2L6 antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using UBE2L6 antibody Catalog No:116531 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IgG (red).

IP Result of anti-UBE2L6 (IP:Catalog No:116531,

UBE2L6 antibody



Catalog Number: 116531

4ug; Detection:Catalog No:116531 1:1000) with Jurkat cells lysate 2000ug.