

# TEAD1 antibody

Catalog Number: 115931

## Product name

TEAD1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TEAD1 Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Positive IP detected in mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 52 kDa

## Preparation

This antibody was obtained by immunization of TEAD1 recombinant protein (Accession Number: BC026959). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:115931 (TEAD1 antibody) at dilution of 1:300



IP Result of anti-TEAD1 (IP: Catalog No:115931, 4ug; Detection: Catalog No:115931 1:300) with mouse skeletal muscle tissue lysate 220ug.