

# Anti-EYA2 Rabbit antibody



Catalog Number: 167931

## Product name

Anti-EYA2 Rabbit antibody

## Specificity

Human

## Antibody description

EYA2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human EYA2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, IP

## Dilutions

WB: 1/2000-1/50000

IP: 1/20

## Validations



Western blot detection of EYA2 in HeLa,A549,U2OS using EYA2 Rabbit mAb(1:1000 diluted)