

NR1H2 antibody



Catalog Number: 113231

Product name

NR1H2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NR1H2 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse liver tissue. Observed molecular weight by Western-blot: 50 kDa, 40 kDa

Preparation

This antibody was obtained by immunization of NR1H2 recombinant protein (Accession Number: NM_007121). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113231(NR1H2 antibody) at dilution of 1:500