

# Human BST2 (His Tag) recombinant protein



Catalog Number: 504231

## General Information

### Gene Name Synonym

HM1.24 antigen; Tetherin

### Protein Construction

A DNA sequence encoding the human BST2 (Q10589) (Asn49-Ser160) was expressed, with an N-terminal polyhistidine tag.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 95 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at -70°C

### Predicted N terminal

His

### Molecular Mass

The recombinant human BST2 comprises 132

amino acids and has a predicted molecular mass of 14.9 kDa. The apparent molecular mass of the protein is approximately 23-26 kDa in SDS-PAGE under reducing conditions due to glycosylation.

### Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human TETHERIN / BST2 / CD317 Protein (His Tag) SDS-PAGE