

IFRD2 antibody



Catalog Number: 111631

Product name

IFRD2 antibody

Specificity

Human; other species not tested.

Antibody description

IFRD2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A549 cells, PC-3 cells. Observed molecular weight by Western-blot: 50 kDa

Preparation

This antibody was obtained by immunization of IFRD2 recombinant protein (Accession Number: NM_006764). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111631(IFRD2 Antibody) at dilution of 1:600