

Rat CLEC4D (Fc Tag) recombinant protein



Catalog Number: 502631

General Information

Gene Name Synonym

C-type lectin superfamily member 8

Protein Construction

A DNA sequence encoding the rat CLEC4D (Q69FH1) (Trp48-Lys218) was expressed with Fc region of human IgG1 at the N-terminus.

Organism

Rat

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu

Molecular Mass

The recombinant rat CLEC4D comprises 431

amino acids and predicts a molecular mass of 48.4 kDa. The apparent molecular mass of the recombinant protein is approximately 53 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Rat CLEC4D / CLECSF8 Protein (Fc Tag) SDS-PAGE